

NEW Physiolaser olympic 500




Physiolaser olympic (with touch screen), laser probe with large touch screen. Available with therapeutic wavelengths up to 500 mW/810 nm cw and 90 W/904 nm pulse 200 nsec..

The new Physiolaser olympic offers more than ever and can be used in all areas of LLLT and biostimulation/biomodulation. It is designed for modern laser therapy, acupuncture + auriculo medicine. The touch screen offers easy operation and fast selection of proven therapeutic biofrequencies:

- NOGIER, high potency up to 18688 Hz
- BAHR 1, high potency up to 38368 Hz
- BAHR 2
- Meridian frequencies
- Chakra frequencies
- All REININGER frequencies
- Free frequency selection and much more

Nogier frequencies and body zones



Many special frequency combinations and therapy programs available. It also offers a sophisticated diagnostic module for point detection and analysis of acupuncture points and organs using the VOLL and NOGIER methods.

Important: Slot for SD memory card to easily upgrade and service the Physiolaser olympic.



Delivery contents:

Main unit, power supply, safety goggles, warning sign, therapy manual, technical manual, safety lock, key, carrying case.

Optional: Probes, applicators

Laser class 3B



NEW Physiolaser olympic 500

Fits in one hand, made for mobile operation with rechargeable battery, versatile and easy to operate.


- Emergency switch
- Safety plug
- Key switch
- Touch screen
- Memory card (update)
- Laser channels/diagnosis
- Laser test



You can use two laser probes simultaneously (single/cluster) and a diagnostic probe in addition (point and organ measurement).

The laser and diagnostic probes offer the highest level of function with display and buttons for one-hand operation (e.g. for RAC).

Display Buttons for one-hand operation



Extensive selection of accessories are available including singleprobes, cluster probes in different configurations, Scanner for larger surfaces and LightNeedle device.




Standard probe with *photon reflection foil*



Brush laser probe for contact therapy



Derma probe with special lens



12 channel LightNeedle device

Please contact us for more information:



RJ-LASER, REIMERS & JANSSEN GmbH
 Medical — Laser — Technology
 Frohnacker 8, 79297 Winden, Germany
 contact@rj-laser.com www.rj-laser.com



Quality management ISO 13485:2003

Distributor